


SSRL  
BL 7-1 MONO

7-7-05



	X	Y	Z
A	-296	-9.948	706.322
B	-300	-9.862	715.574
C	-16.736	-9.955	706.322
D	-16.749	-9.861	715.581

SSRL  
BL 7-1 MONO  
LAB DATA

7-7-05  
J.M., L.G.

(X)

CRYSTAL

W/S) 11.509 READ ON (+) SIDE  
.110 (.118 X COS 21.88')

11.619

(-13.001) VALUE

+8.618 LOS

A) 7.914  
1.000

8.914

8.618

-2.96

B) 7.918  
1.000

8.918

8.618

-.300

C) 16.940 (A-C) D) 24.367  
1.500 T/B 1.000

16.440

-.296

-16.736

25.367

8.618

-16.749

(Z)

708.936 W/S END  
11.806 OF CRYSTAL

697.130 LOS

A) 8.192  
1.000

9.192

697.130

706.322

B) 17.444  
1.000

18.444

697.130

715.574

C) 8.192  
1.000

9.192

697.130

706.322

D) 17.451  
1.000

18.451

697.130

715.581

SSRL  
 B/L 7-1 MONO  
 LAB DATA check  
 Jm

7-7-05  
 J.M., L.G.

(Y)

CRYSTAL

u/s) +.024 READ	D/S) +.085
.118	.118
+ .094 FROM $\epsilon$	+ .033 FROM $\epsilon$
2.696 VAC. @ $\epsilon$	2.755
2.790	2.788
M = 2.789 Hz	




A) 20.545	B) 20.459	C) 20.552	D) 20.458
1.250	1.250	1.250	1.250
21.795	21.709	21.802	21.708
11.847	11.847	11.847	11.847
-9.948	-9.862	-9.955	-9.861

GUN1	GUN2
9.578 REF.	18.636 GUN2
2.789 Hz	6.789 $\epsilon$
6.789 $\epsilon$	+11.847 Hz


BL7-1 Mono

Alignment

Plan View




Elevation Along Optical Path (5.4 mrad)



### Rotation Axis Alignment

1. Optical Face of Crystal must be parallel to Rotation Axis such that the crystal face will sweep a plane that is orthogonal to the rotation axis. This can be achieved by aligning the crystal pitch and roll parallel to the Huber Plate

set screws for planer alignment




Huber Plate

A. Roll w.r.t. Huber Plate: Best accuracy - limited to measuring geometry of crystal or Flexure Holder  
 $\pm 2 \text{ mrad}$

B. Pitch, measured over 7 inches  $\pm .005^\circ / \text{in} = \pm 0.7 \text{ mrad}$   
(w.r.t. Huber)

C.  $\psi$  w.r.t. Huber Plate set to  $\pm 0.1^\circ$   
measured  $\pm .002^\circ$

2. Crystal Yaw (Theta) to be set in known angle. By setting indicated edge of glid cap parallel to wiggler  $\Phi$  in X, the crystal should be at 21.88 deg w.r.t. wiggler  $\Phi$



$$\pm 0.01/5^\circ = \pm 2 \text{ mrad}$$

Make this edge parallel to wiggler  $\Phi$  in X. ( $X = -2.496$ )

Edge of Huber Plate is parallel to wiggler  $\Phi$  in X.

### 3. Beam Line Alignment

Once the above conditions are set, the whole assembly can be set into installation coordinates. This involves tilting the rotation axis  $9 \text{ mrad}$  as indicated on page 1.


Y & X become critical for alignment.

Pitch  $9 \text{ mrad} \pm 0.3 \text{ mrad}$

X @  $1/3$  end of crystal  $\pm .005$

Y @  $1/3$  end of crystal  $\pm .002$

Target B




BL 7-1 Mono

Alignment

3. Beam Alignment - Continue

Target A & B @  
E of crystal



2.654 D


Target A  $x = -3.001''$   $y = 2.696''$

$z = 708.936''$

Target B  $x = -5.657''$   $y = 2.755''$

$z = 715.547''$

$\Delta y = 1.059$



SSRL  
BL7-1 MONO  
LAB DATA

7-7-05

(X)


check

GLIOLOP

11.100 READ ON GLIOLOP w/s & p/s PARALLEL TO HUBER PLATE  
- 2.496 VAL.  
+ 8.604 LOS

CRYSTAL

u/s  
11.509 READ w/ (+) SIDE OF CRYSTAL  
.110 (118X4521.BB)  
11.619 to (-) SIDE  
8.604  
3.015



+ 8.604 LOS  
10.000 REF  
18.604 = 4

D/S

- 5.657  
.110  
5.547  
18.604  
24.151

3.015 + 0.14  
5.557 D/S + 0.10

$\lambda = 21.865^{\circ}$

0 SIR  
+ 0.10

$D = 2.652$

A-C  
16.940